




DISINFECTANT
SOLUTIONS

Triple Quick

Fresh & Clean
Disinfecting Cleaner



- Broad-spectrum power kills a wide range of pathogens
- Versatile formula cleans and disinfects for one-step use
- Ready-to-use formula requires no diluting or mixing

State[®]
Since 1911

stateindustrial.com

Care For Work Environments[®]

Triple Quick

Fresh & Clean
Disinfecting Cleaner

PRODUCT OVERVIEW :

CLEANER, DISINFECTANT AND DEODORIZER

This ready-to-use product prevents cross contamination on facility surfaces. Contains biodegradable surfactants, wetting agents and powerful sequestering agents to provide maximum cleaning power.

READY-TO-USE

Handy, ready-mixed, quart size eliminates the need for measuring and mixing.

WHERE TO USE:

Effective on all hard, non-porous, inanimate surfaces such as floors, walls, metal surfaces, stainless steel surfaces, porcelain, glazed ceramic tile, plastic surfaces, bathrooms, shower stalls, bathtubs and cabinets.

TUBERCULOCIDAL

Triple Quick is effective against Mycobacterium tuberculosis (BCg), which can cause Tuberculosis. Triple Quick kills HIV-1 (AIDS virus) in hospitals and other health care facilities.

EFFECTIVE AGAINST CANINE PARVOVIRUS

Use on inanimate, hard, non-porous surfaces such as kennel walls and floors, and veterinary examination tables and counters. Triple Quick kills the Canine Parvovirus, which is a potentially fatal disease for dogs and other members of the canine family.

BACTERICIDE

Effective against Staphylococcus aureus, Salmonella choleraesuis, and Pseudomonas aeruginosa. Triple Quick meets the strict requirements for hospital use.

FUNGICIDE

Kills pathogenic fungi (such as T. mentagrophytes, the fungus that causes Athlete's foot) on hard surfaces found in bathrooms, shower stalls, locker rooms and other clean, non-porous, hard surfaces commonly contacted by bare feet.

MILDEWSTAT

Controls mold and mildew effectively. Apply Triple Quick at weekly intervals or when mildew growth appears.

DIRECTIONS FOR USE :

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body.

For plastic and painted surfaces, spot test on an inconspicuous area before use.

TECHNICAL DATA :

APPEARANCE: clear, pale yellow liquid

SCENT: Fresh & Clean

WEIGHT PER GALLON: 8.4 lbs/gal

pH: approximately 12.0

DISINFECTANT TYPE: quaternary ammonium compound

SHELF STABILITY: two years

REGISTRATIONS: DIN 02428377

AUTHORIZATIONS: Meets the requirements for use in federally inspected meat and poultry plants.

For more product information and SDS, visit stateindustrial.com.

PACKAGING

125843

1 QT Bottle/CS12

CONTACT US

 stateindustrial.com

 support@stateindustrial.com

 1.866.747.2229

 5915 Landerbrook Drive, Mayfield Heights, OH 44124

