

DISINFECTANT *Solutions*



- Hospital-use and broad-spectrum disinfectant when used as directed
- Unique locking spray valve allows for total-release of product for an adjunct to normal disinfection
- Pleasant citrus fragrance quickly eliminates odors from many sources as a total release product



CARE FOR WORK ENVIRONMENTS®



Total Release

VERSATILE DISINFECTANT

Discover TR is a disinfectant, antimicrobial agent and two-way deodorizer for use on pre-cleaned, hard, non-porous and non-food contact surfaces. Discover TR is a virucidal, tuberculocidal, fungicidal, bactericidal, Pseudomonacidal, and Staphylocidal. It is designed for use in/on industrial, institutional, commercial, medical, and residential facilities, equipment, and vehicles which can harbor disease causing microorganisms. Fogging is an adjunct or supplement to normal cleaning and disinfection procedures and practices.

VIRUCIDAL

Kills Human Immunodeficiency Virus Type 1 (HIV-1; a retrovirus that causes AIDS), Herpes Simplex Virus Type 1 and Herpes Simplex Virus Type 2, Influenza A2/Hong Kong Virus, Adenovirus Type 5, Vaccinia Virus, Infectious Canine Hepatitis Virus, Canine Distemper Virus, Feline Pneumonitis Virus, and Avian Influenza A Virus (H5N1). Follow label instructions for disinfection.

BACTERICIDAL

Kills Escherichia coli; Pseudomonas aeruginosa; Salmonella enterica; Salmonella schottmuelleri (paratyphoid B); Shigella dysenteriae; Streptococcus pyogenes; Staphylococcus aureus; and Methicillin-resistant Staphylococcus aureus (MRSA).

TUBERCULOCIDAL

Kills Mycobacterium tuberculosis bovis BCG (Tubercle bacilli). Follow label instructions for disinfection.

FUNGICIDAL

Kills Trichophyton mentagrophytes (athlete's foot fungi; formerly known as T. interdigitale) in bathrooms, shower areas, and locker rooms. Follow label instructions for disinfection.

CONTROLS MOLD/MILDEW

Completely inhibits the growth of Aspergillus niger and Penicillium glaucum fungi on fabric.

SANITIZES

Controls Staphylococcus aureus and Klebsiella pneumoniae on porous fabric and leather surfaces, such as on shoes, athletic mats, athletic shoes and other athletic equipment.

DEODORIZES

Pleasant fragrance, and disinfectant quickly eliminates odor-causing bacteria.

TOTAL-RELEASE FOGGER

Discover TR is packaged in a unique, total-release can that creates an enveloping fog for rooms and spaces greater than 5 ft. X 5 ft. Fogging is an adjunct or supplement to normal cleaning and disinfection procedures and processes.

PACKAGING

121222-12 6 OZ Aerosol/CS12
121222-24 6 OZ Aerosol/CS24

FOR USE IN

Industrial, institutional, commercial, medical, and residential facilities, equipment, and vehicles, including hospitals ambulances, nursing homes, medical and dental offices and clinics, veterinary offices and animal areas, schools, laboratories, hotels and motels, restaurants and cafeterias, factories and manufacturing plants, and apartments buildings and homes.

FOR USE ON

Most glass, glazed porcelain and ceramic, plastic, rubber, metal, finished wood, sealed leather, sealed grout, and similar nonporous surfaces.

DIRECTIONS FOR FOGGING ROOMS GREATER THAN 5 FT X 5 FT

All surfaces must be cleaned and disinfected in accordance with label directions prior to fogging.

Discover TR can only be used in temperatures greater than 50°F. Do not use more than one fogger per room. Do not use in small, enclosed spaces such as closets, cabinets or under counters or tables. Do not use in a room 5 ft. x 5 ft. or smaller. Turn off ALL ignition sources such as pilot lights (shut off gas valves), or open flames or running electrical appliances that cycle off and on (i.e. refrigerators, thermostats, etc.). Call your gas utility or management company if you need assistance with your pilot lights. Remove exposed food, utensils, and dishes. Cover immobile food handling equipment. Open cabinets and doors to area to be treated. Shut off fans and air conditioners. Seal doors and windows.

Use one unit for each 6,000 cubic ft. of unobstructed area. Do not remain in the area during treatment. Ventilate thoroughly before re-entry. Do not use this unit in a cabinet or under a counter or table. Do not use this in an area less than 100 cubic ft. Do not use more than one unit per average size room. Point valve opening away from face and eyes when releasing.

TO OPERATE VALVE: Place several layers of newspaper or pad on stand or table in the center of room. To lock fogger valve in open position for automatic discharge, press valve button all the way down hooking in catch. Then place fogger on the stand or table. Leave room or building and keep room/building closed for two hours before airing out. Open all doors and windows and allow to air for five hours. Do not attempt to unlock valve once fogging has started.

TECHNICAL DATA

Appearance: Clear, colorless aerosol spray

Fragrance: Citrus

Disinfectant Type: Phenolic

Active Ingredients: 2-phenylphenol 0.209% -- 4-Tert-Amylphenol 0.05% -- Ethanol 53.00%

EPA Reg. No. 44446-67-70799

EPA Est. No. 44446-TX-1



State Industrial Products
5915 Landerbrook Drive
Mayfield Heights, OH 44124
To Order Call: 1-866-747-2229
www.statechemical.com

Puerto Rico
Royal Industrial Park, Bldg "M"
Local #5, Carr 869, km 1.5 Palmas
Cataño, P.R. 00962
To Order Call: 787-275-3185

