



DISINFECTANT SOLUTIONS

Quick DEFENSE™ NDC



- Super concentrated, pH-neutral, broad-spectrum quaternary disinfectant
- Highly effective against SARS-CoV-2, HIV-1 virus, HCV, Herpes, MRSA, and more
- Cleans and disinfects hard, non-porous surfaces



stateindustrial.com

Care For Work Environments®



Quick DEFENSE™ NDC™

PRODUCT OVERVIEW :

FAST-ACTING SOLUTION

Quick Defense™ NDC™ is a super concentrated, pH-neutral, broad-spectrum quaternary disinfectant that is specially formulated to clean and disinfect virtually any hard, non-porous surface. Its unique odor elimination kills a variety of odor-causing microorganisms. Its unique formula works twice as fast on most pathogens as other products the market, making Quick Defense™ NDC™ the ideal solution for healthcare applications and restrooms.

MARKETS

Quick Defense™ NDC™ is the perfect disinfectant solution for a variety of markets. Utilize Quick Defense™ NDC™ in healthcare facilities or institutions, correctional facilities, hotels/motels, schools, animal care facilities, nursing homes, daycare centers, restaurants, funeral homes and recreational facilities.

KILL CLAIMS

Quick Defense™ NDC™ is a super concentrated formula that is highly effective against dangerous bacteria and viruses such as:

- SARS-CoV-2
- Salmonella
- Herpes Simplex Virus Types 1 & 2
- MRSA
- HBV (Hepatitis B Virus)
- Staphylococcus aureus
- E.coli
- HCV (Hepatitis C Virus)
- Streptococcus pyogenes
- Listeria monocytogenes
- HIV-1 (AIDS Virus)

For use against SARS-CoV-2, the cause of COVID-19. Quick Defense™ NDC™ has a 2-minute dwell time.

COMMON APPLICATIONS

The quaternary, pH-neutral formula of Quick Defense™ NDC™ will not harm or dull high-gloss floor finishes. It is ideal for use on:

- Floors
- Metals
- Glazed ceramic tile
- Showers
- Toilets and urinals
- Mirrors
- Stainless steel
- Cabinets
- Glass
- Glazed porcelain
- Plastics
- Stalls

DIRECTIONS FOR USE :

For general disinfection, dilute at a ratio of 1 oz. per 1 gallon of water. See product label for correct dwell times.

To comply with federal law, use this product in accordance with its applications, directions and procedures on the label.

Please review the label and SDS on stateindustrial.com for all product directions, precautions, and first aid information.

TECHNICAL DATA :

APPEARANCE: colorless liquid

FRAGRANCE: orange

DISINFECTANT TYPES: quaternary compounds

pH: 6.0-8.0


PACKAGING

127889	20 GL Drum
127888	5 GL Pail
127887	1 GL Bottle/CS4
127885	3 L Bottle/CS4

CONTACT US

 stateindustrial.com

 support@stateindustrial.com

 1.866.747.2229

 5915 Landerbrook Drive, Mayfield Heights, OH 44124

